

Meas. No. 7511 SOP 45 Date 07.09.2009 16:24:52 Operator: FRITSCHLAN\Benes

Material: Dent. compos. filler (fine)

P7pl ZrO2 80ml 2mm 4h

Description: 2min ultrasonic + water + 0,1%Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 9,9 %

Pump 65,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090188**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
1,137	5,0
1,369	10,0
1,540	15,0
1,690	20,0
1,831	25,0
1,969	30,0
2,107	35,0
2,260	40,0
2,423	45,0
2,601	50,0
2,785	55,0
3,017	60,0
3,287	65,0
3,632	70,0
4,162	75,0
5,035	80,0
6,355	85,0
7,774	90,0
9,438	95,0
12,086	99,0

01-20µm

Obere Kornklasse [µm]	Q3(x) [%]
0,100	0,0
0,150	0,0
0,200	0,0
0,300	0,3
0,400	1,0
0,500	1,7
0,600	2,0
0,700	2,2
0,800	2,3
1,000	3,3
1,200	6,1
1,500	13,7
2,000	31,1
2,500	47,2
3,000	59,6
4,000	73,9
5,000	79,8
6,000	83,7
8,000	90,8
10,000	96,1
15,000	99,9
20,000	100,0

