

Meas. No. 7532 SOP 45 Date 09.09.2009 14:47:27 Operator: FRITSCHLAN\Benes

Material: Dent. compos. filler (fine)

P7pl ZrO2 80ml 2mm 8h

Description: 2min ultrasonic + water + 0,1%Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 10,9 %

Pump 65,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090188**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,708	5,0
0,954	10,0
1,169	15,0
1,378	20,0
1,596	25,0
1,829	30,0
2,082	35,0
2,382	40,0
2,720	45,0
3,112	50,0
3,580	55,0
4,175	60,0
4,904	65,0
5,785	70,0
6,858	75,0
8,104	80,0
9,593	85,0
11,520	90,0
14,341	95,0
19,506	99,0

01-30µm

Obere Kornklasse [µm]	Q3(x) [%]
0,120	0,6
0,200	1,4
0,300	1,6
0,500	2,3
1,000	11,0
1,500	22,9
2,000	33,4
3,000	48,6
4,000	58,7
5,000	65,6
6,000	71,0
7,000	75,6
8,000	79,6
9,000	83,1
10,000	86,1
15,000	95,8
20,000	99,2
25,000	99,9
30,000	100,0

