


% <	μm [ø]	CV [%]	M89226	M89227
10,0	11,65	1,35	11,81	11,5
50,0	29,20	0,34	29,3	29,1
90,0	154,38	3,61	159,96	148,81

< μm	% [ø]	CV [%]	M89226	M89227
0,10	0,5	59,99	0,2	0,8
1,00	2,3	14,09	2	2,6
5,00	3,8	6,55	3,5	4
10,00	7,7	4,17	7,4	8
15,00	16,2	0,24	16,1	16,2
20,00	28,0	0,54	28,1	27,8
30,00	51,7	0,51	51,5	52
50,00	75,9	1,20	75	76,8
80,00	79,7	1,11	78,8	80,6
100,00	81,1	1,15	80,1	82
150,00	89,3	1,05	88,3	90,2
200,00	96,0	0,67	95,4	96,6
300,00	99,9	0,05	99,8	99,9
400,00	100,0	0,00	100	100
500,00	100,0	0,00	100	100



HDPE "N+S"

Mess.-Nr. 89226 Datum Freitag, 20 Februar 2015 11:39 Bearbeitungs- Nr.:

Material HDPE "N+S"

Bemerkung Trockenmessung, Druck: 2 bar, Rinne: 2,5 mm


Berechnung Fraunhofer broad TradeOff 1,000000e+003

Scans 30 Scans Grobbereich 30 Messkanäle 102 Median 29,27 [µm]

Messbereich 0,08 [µm] - 2000,00 [µm]

%	<µm
5	7,64
10	11,81
15	14,47
20	16,72
30	20,72
40	24,78
50	29,30
60	34,69
70	42,64
80	98,50
90	159,96
95	196,68
97	218,27
99	255,14

%	<µm
0,2	0,1
2,0	1,0
3,5	5,0
7,4	10,0
16,1	15,0
28,1	20,0
51,5	30,0
75,0	50,0
78,8	80,0
80,1	100,0
88,3	150,0
95,4	200,0
99,8	300,0
100,0	400,0
100,0	500,0




89226 dQ3(x) — 89226 Q3(x)

Mess.-Nr. 89227 **Datum** Freitag, 20 Februar 2015 11:40 **Bearbeitungs- Nr.:**
Material HDPE "N+S"
Bemerkung Trockenmessung, Druck: 2 bar, Rinne: 2,5 mm

Berechnung Fraunhofer broad **TradeOff** 1,000000e+003
Scans 30 **Scans** 30 **Messkanäle** 102 **Median** 29,07 [µm]
Feinbereich **Grobbereich**
Messbereich 0,08 [µm] - 2000,00 [µm]

%	<µm
5	6,77
10	11,50
15	14,39
20	16,78
30	20,84
40	24,81
50	29,10
60	34,07
70	41,05
80	63,45
90	148,81
95	183,33
97	203,86
99	240,23

%	<µm
0,8	0,1
2,6	1,0
4,0	5,0
8,0	10,0
16,2	15,0
27,8	20,0
52,0	30,0
76,8	50,0
80,6	80,0
82,0	100,0
90,2	150,0
96,6	200,0
99,9	300,0
100,0	400,0
100,0	500,0



 89227 dQ3(x)
  89227 Q3(x)